PANSERIN™ 604

PANSERIN™ 604 is a complete medium ready for use for the cultivation of transfected and non-transfected CHO-cells (Chinese Hamster Ovary). These cells are frequently used for the production of recombinant proteins. The cells normally grow adherently (e.g. in roller cultures), however with suitable media and with the optimal media they can also be kept in suspension culture (continuous perfusion or batch systems).

PAN-Biotech has developed an optimal version for every application:

For the suspension culture:

- **PANSERIN™ 604ST serumfree** – optimal cell growth, quick adaption
- **PANSERIN™ 604S without protein** or peptide components of **animal** or **human origin**
- **PANSERIN™ 604SPF** **chemically defined, proteinfree**

The following graph shows the growth curves of the PANSERIN™ 604 versions.

The CHO-cells were cultivated from an adherent culture containing serum in PANSERIN™ 604 without a preceding adaption. The seed amounted to 50 000 cells/ml. The maximum generation times were 15,3 to 18,5 hours.

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
PAN[™]
BIOTECH GmbH

Composition:**PANSERIN™ 604 serumfree:**

A basal medium optimized for suspension cultures, which had been designed for a high cell density, was supplemented with hormones, lipids and special growth factors for an optimal growth of CHO-cells. It does not contain any undefined lysates or plant hydrolysates. It allows a very quick growth of the cells without long adaption phases.

PANSERIN™ 604S without animal or human components

Basal medium which has been completed with exclusively recombinant proteins. It does not contain any undefined lysates or plant hydrolysates. Especially suitable for the production of proteins which are applied in sensitive fields. It excels by a quick cell growth without long adaption phases.

PANSERIN™ 604SPF proteinfree, chemically defined:

This medium is completely proteinfree and chemically defined. After a short adaption phase excellent growth rates can be reached with this medium and these growth rates are only slightly under the culture media containing protein and serum. Especially suitable for the production of recombinant proteins for the application in sensitive fields.

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
**CHO-Zellen in PANSERIN™ 604ST, serumfree**

CHO-Cells were cultivated directly without adaption phase from an adherent culture containing serum in a suspension culture with PANSERIN™ 604 (shaking culture).

Because of an optimal and well-balanced composition of the medium an aggregate formation of the cells, as it can be seen with many serumfree media, could be largely avoided (see illustration).

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PAN™
BIOTECH GmbH

For the adherent culture:

PANSERIN™ 604 supports the growth of adherent CHO-cells in an optimal way. With attachment factors the cells can grow on quickly and this allows a high cell output within a short time.

**CHO-Cells in PANSERIN™ 604**

CHO-Cells were transferred from a culture containing serum directly into PANSERIN™ 604.

Composition:

A basal medium optimized for CHO-cells was supplemented with albumin, insulin, transferrin, attachment factors, growth factors, lipoproteins and trace elements. It does not contain any undefined hydrolysates or lysates.

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
CHO-Cells were cultivated from a culture containing serum in PANSERIN™ 604 without adaption. As a comparison CHO-cells in DMEM with 4,5 g/l glucose, L-glutamine and pyruvate, supplemented with 10 % FBS.

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